

This is a digital copy of a book that was preserved for generations on library shelves before it was carefully scanned by Google as part of a project to make the world's books discoverable online.

It has survived long enough for the copyright to expire and the book to enter the public domain. A public domain book is one that was never subject to copyright or whose legal copyright term has expired. Whether a book is in the public domain may vary country to country. Public domain books are our gateways to the past, representing a wealth of history, culture and knowledge that's often difficult to discover.

Marks, notations and other marginalia present in the original volume will appear in this file - a reminder of this book's long journey from the publisher to a library and finally to you.

### Usage guidelines

Google is proud to partner with libraries to digitize public domain materials and make them widely accessible. Public domain books belong to the public and we are merely their custodians. Nevertheless, this work is expensive, so in order to keep providing this resource, we have taken steps to prevent abuse by commercial parties, including placing technical restrictions on automated querying.

We also ask that you:

- + *Make non-commercial use of the files* We designed Google Book Search for use by individuals, and we request that you use these files for personal, non-commercial purposes.
- + Refrain from automated querying Do not send automated queries of any sort to Google's system: If you are conducting research on machine translation, optical character recognition or other areas where access to a large amount of text is helpful, please contact us. We encourage the use of public domain materials for these purposes and may be able to help.
- + *Maintain attribution* The Google "watermark" you see on each file is essential for informing people about this project and helping them find additional materials through Google Book Search. Please do not remove it.
- + *Keep it legal* Whatever your use, remember that you are responsible for ensuring that what you are doing is legal. Do not assume that just because we believe a book is in the public domain for users in the United States, that the work is also in the public domain for users in other countries. Whether a book is still in copyright varies from country to country, and we can't offer guidance on whether any specific use of any specific book is allowed. Please do not assume that a book's appearance in Google Book Search means it can be used in any manner anywhere in the world. Copyright infringement liability can be quite severe.

#### **About Google Book Search**

Google's mission is to organize the world's information and to make it universally accessible and useful. Google Book Search helps readers discover the world's books while helping authors and publishers reach new audiences. You can search through the full text of this book on the web at http://books.google.com/



BRIEF COURSE IN

# METALLURGICAL ANALYSIS

BY

HENRY ZIEGEL

# Library

of the

University of Wisconsin

٠		
•		
		4

		•	
-			

METALLURGICAL ANALYSIS.

#### Published by

# The Chemical Publishing Co.

Easton, Penna.

Publishers of Scientific Books

Engineering Chemistry

Portland Cement

Agricultural Chemistry

**Qualitative Analysis** 

Household Chemistry

Metallurgy, Etc.

# BRIEF COURSE IN

# Metallurgical Analysis

BY

Henry Ziegel

TO PROFESSOR G. P. BAXTER, Teacher, Investigator, Friend

EASTON, PA.: THE CHEMICAL PUBLISHING CO. 1915

LONDON, ENGLAND:
WILLIAMS & NORGATE
14 MENRIETTA STREET, COVENT GARDEN, W. C.

COPYRIGHT, 1915, BY EDWARD HART.

209253 MAR -3 1917 MKT .ZG

# CONTENTS.

	PAGE
Introduction	V
Iron Ores, Iron and Steel·····	1
Silica in Ores	I
Silicon in Iron and Steel	2
Phosphorus Determinations, as Yellow Precipitate · · · · · · · · · · · · · · · · · · ·	2
Sulfur Determinations, Inclusion, Occlusion	6
Carbon by Direct Combustion	IO
Manganese by the Bismuthate Method	12
Conditions for the Titration of Iron by Permanganate	14
Iron Ores, Iron and Steel, Lesser Constituents	20
Nickel Determinations as Dimethylglyoxime	20
Copper in Steel, Electrolyis	22
Tungsten	22
Chromium in Ores and Steel·····	24
Vanadium in Ores and Steel	24
Chromium and Vanadium in Steel. Usual Method	26
Molybdenum	26
Titanium in Iron	28
Aluminum	30
Limestones, Slags and Clays	32
Limestone, Determination of Insoluble, Al <sub>2</sub> O <sub>3</sub> , Fe <sub>2</sub> O <sub>3</sub> , CaO, MgO	32
CaO and MgO in same solution	34
Mg electrometrically	34
Slags, Determination of SiO <sub>2</sub> , Cu, Fe, CaO	34
Na <sub>2</sub> O and K <sub>2</sub> O in Clays	36
Non-Ferrous Alloys	38
Brass, Determination of Copper, Lead, Zinc, Tin, Iron, Nickel, discussion of methods	38
Bronze, SnO <sub>2</sub> purification	42
New Rapid Analysis	46
Phosphorus and Manganese in	46
Other similar alloys	48
Bearing Metals Analysis of Lead, Antimony, Tin, Copper, Zinc	48
Non-Ferrous Ores, Matte	54
Matte	54
Ores	58
Copper in Ores	58
Zinc in Ores	60
Lead in Ores	62
Al <sub>2</sub> O <sub>2</sub> in Bauxite	66
Pibliography	60

		:
·		
		!

# INTRODUCTION.

This book is written for the use of upperclassmen, trained in elementary qualitative and quantitative analysis. The work outlined in the procedures should require one-fourth to one-fifth the students time, some ten to twelve hours a week—during two successive semesters. Since few calculations are implied, the student should have or have had accompanying his previous course, a series of problems in chemical arithmetic.

The author's object then, has not been to describe all the methods for determining the composition of metallurgical products, rocks and ores, but to describe a series of standard analyses in daily use in various analytical laboratories. The procedures have been so chosen so as to avoid the same determinations as much as possible. Procedures for similar substances have been mentioned in the notes or referred to in the bibliography. Since such a selection is always open to comment, the author welcomes criticism and suggestion.

The analyses of gases, fuels and water, form a sub-division which may well be treated in a separate course; but several articles on sampling<sup>1</sup> should be read by the student.

By binding the laboratory sheets with the notes, neatness and original data are continually suggested to the student.

Finally, due credit must be given to the authors of the larger works on the subject, vis.:

Blair: "The Chemical Analysis of Iron." 7th ed.

Johnson: "The Chemical Analysis of Special Steels." 2nd ed.

Low: "Technical Methods of Ore Analysis."

Lord and Demorest: "Metallurgical Analysis."

Price and Meade: "Technical Analysis of Brass."

Smith: "Electro-Analysis" 5th ed.

Classen: "Quantitative Analysis by Electrolysis," Tr. Hall.

Betts: "Lead Refining by Electrolysis."

"Methods of Analysis of Carbon Free Metals," Monograph by Goldschmidt Co.

"Treadwell," 3rd ed., and "Fresenius," 6th ed., on the entire subject of analysis.

The Journals and Abstracts of The American Chemical Society.

The publications of the American Electrochemical Society.

The publications of the Bureau of Standards, Washington.

The publications of the Bureau of Mines, Canada.

The reader's attention is called to the procedures and notes on the analyses including the use of the hydrogen and calomel electrode<sup>2</sup> (work of Sand, Schoch, 'Hillebrand and Forbes and Bartlett). To be sure these procedures require special equipment and some skill, yet the results obtained with them are so satisfactory and lead to such interesting future application that the student should study these procedures in detail.

# IRON ORES, IRON AND STEEL, CHIEF CONSTITUENTS.

# SILICA IN ORES.'

There are two general methods for the determination of silica in ores; by solution and by fusion. Several general details of manipulation are emphasized once for all in the first few procedures.

Procedure No. 1.—Dissolve I gram of ore in 15 cubic centimeters HCl specific gravity 1.2 in a small beaker provided with a cover glass; warming until the ore is decomposed. Oxidize any ferrous iron with a few drops of HNO<sub>3</sub>, continue heating, wash off the cover glass, and evaporate until the solution is quite syrupy. Evaporate to dryness and bake at 110° C. for 15 minutes. Cool, redissolve in HCl, evaporate and rebake. Redissolve in 30 cubic centimeters HCl 1-3, and filter through a small ashless paper using a policeman of ½-inch glass rodding tipped by a small piece of well washed rubber tubing. It is well to remember that many applications of small quantities of water will wash out more precipitate than few applications of less small quantities.

Wash the filter paper free from iron by using warm wash water containing a little HCl (5 to 10 cubic centimeters per liter), remove paper from funnel with tip of finger nail or drawn out end of a policeman, and place in a weighed 20-gram platinum crucible. Char, ignite, cool and weigh. Add two drops H<sub>2</sub>SO<sub>4</sub> specific gravity 1.84, and 7 cubic centimeters H<sub>2</sub>F<sub>2</sub>. Ignite over a ring burner,\* cool and reweigh. Difference between weighings is SiO<sub>2</sub>. Wait at least 30 minutes after putting crucible into desiccator before removing to weigh.

Procedure No. 1, part 2.—Place I gram of ore with 7 to 8 grams sodium carbonate in a 30-gram platinum crucible, stir the mixture thoroughly with a thin glass rod, wipe off latter with a piece of filter paper, and add the paper to the charge. Heat to quiet fusion, chill the crucible in water, and dissolve

<sup>\*</sup>The cylinder described by Hildebrand, Bull. U. S. G. S., No. 422, p. 31, may be used instead.

the melt in 40 cubic centimeters HCl 1-3. Remove the crucible and cover from the HCl solution with a glass rod, add a few drops HNO<sub>3</sub>, and proceed as above, oxidizing with HNO<sub>3</sub> ff.

## SILICON IN IRON AND STEEL.

Procedure No. 2.—Dissolve from 2 to 5 grams of steel in 20-40 cubic centimeters nitric acid 1-1 in a platinum dish covered with a cover glass, wash off cover glass into dish and evaporate contents to dryness on the air bath. Increase the heat to decompose ferric nitrate, and bake for ten minutes. Redissolve in 10 cubic centimeters concentrated HCl slowly and again evaporate to dryness and bake. Redissolve in HCl and add 100 cubic centimeters water, decant, filter and purify as before. Factor ½ of 0.4693 for 5 grams.

Or silicic acid in dilute acid solution may be completely dehydrated by 20 cubic centimeters of H<sub>2</sub>SO<sub>4</sub> 1-1 upon boiling to fumes of SO<sub>3</sub> cooling and diluting. Filter hot and determine the silica as before.

Twenty-five cubic centimeters of the following solvent may be added all at once. Water, 50 cubic centimeters; concentrated H<sub>2</sub>SO<sub>4</sub>, 40 cubic centimeters; concentrated HNO<sub>3</sub>, 80 cubic centimeters; concentrated HCl, 60 cubic centimeters. Ferro silicons require fusion with Na<sub>2</sub>CO<sub>3</sub> and Na<sub>2</sub>O<sub>2</sub>, preferably in a nickel crucible.

Silica in slag is determined by moistening 5 grams in a casserole with a few drops of water stirring in 5 cubic centimeters HCl and warming until the slag gelatinizes. Two evaporations from HCl solution are followed by the usual fusion and purification.

## PHOSPHORUS.

Phosphorus is determined as the yellow precipitate, (NH<sub>4</sub>)<sub>3</sub> PO<sub>4</sub>. 12MO<sub>3</sub> containing 1.63 per cent. P. The precipitate may be weighed as such, or dissolved in a pipetteful of standardized NaOH, and the excess alkali titrated with standard HNO<sub>3</sub>, using phenolphthalein as an indicator.

The following conditions must be observed for quantitative precipitation. Absence of HCl and of organic matter.

Presence of ammonia as nitrate or sulphate.

Analysis of	Sample No.		
		-	
,			

Presence of 30-60 molecules of sulphuric or better nitric acid per molecule  $P_aO_a$ .

Phosphorus should all be present as tribasic phosphoric acid.

Temperature should never exceed 80° C., and should be around 40° C. (At elevated temperatures MO<sub>8</sub>, which cannot be removed, will precipitate, so also may compounds of iron and ammonia with MO<sub>8</sub>.) A convenient procedure is to neutralize the paste of Fe (OH)<sub>8</sub> with dilute nitric acid, and after adding the precipitating solution, merely set aside to cool and coagulate for about an hour.

Procedure For Iron Ores, No. 3.—Weigh 2, 1.63, or a multiple grams of ore into 1/4-liter Erlenmeyer flask and dissolve in 15-40 cubic centimeters concentrated HCl. Filter and wash; ignite residue. Fuse latter with 2 or 3 grams Na<sub>2</sub>CO<sub>3</sub>, dissolve fusion in HCl, evaporate both solutions down to dryness and bake, preferably repeating. Cool, dissolve in as little HCl as possible, combine the two solutions and filter off SiO<sub>3</sub>.

To the filtrate in an Erlenmeyer flask add 30 cubic centimeters concentrated HNO<sub>3</sub> free from yellow fumes and evaporate to syrupy; repeat so as to rid of last traces of HCl.

(Iron and steel should be dissolved in dilute HNO<sub>3</sub>, the silica dehydrated (as Procedure No. 1); the resulting solution in HCl treated with HNO<sub>3</sub>, and the following common procedure followed.)

Add NH<sub>4</sub>OH to filtrate in an Erlenmeyer flask to form a thick paste, and dissolve this paste in HNO<sub>3</sub> I-I, giving somewhat less than 200 cubic centimeters of amber colored solution at about the proper temperature to add 60 cubic centimeters of the precipitant.\* Shake the solution for five minutes, remove stopper and let stand for I hour. Filter into a weighed gooch crucible, wash precipitate thoroughly with IO per cent. NH<sub>4</sub>NO<sub>3</sub> solution, twice with water and twice with alcohol, dry for 40 minutes at IIO° C., cool and weigh.

<sup>\*</sup> Molybdate solution. Dissolve 10 grams of molybdic anhydrid in 50 cubic centimeters H<sub>2</sub>O containing 8 cubic centimeters NH<sub>4</sub>OH. Filter through asbestos and pour it into a mixture of 40 cubic centimeters HNO<sub>2</sub> and 60 cubic centimeters H<sub>2</sub>O. Stir while mixing, and keep cool, add one drop of microcosmic salt solution, and let settle for 24 hours before decanting for use.

Analysis of Sam	ple No.			
		· · · · · · · · · · · · · · · · · · ·		
	<del></del>			
	·			
<u> </u>				
· · · · · · · · · · · · · · · · · · ·				
				!
		·		
			•	

•

Instead of weighing the yellow precipitate it may be dissolved in a pipette full of NaOH and the excess NaOH titrated with HNO<sub>3</sub> using I drop of a I per cent. alcoholic solution of phenolphthalein as an indicator. In this case, an ore or a steel of known phosphorus content should be run at the same time, so as to determine the phosphorus equivalent of the standard acid and alkali.

Standard HNO<sub>3</sub>, 20 cubic centimeters concentrated HNO<sub>3</sub> per liter of water. Standard NaOH, saturate 30 cubic centimeters H<sub>2</sub>O with NaOH and use upper 3/3 in a liter of water. These solutions should be made equivalent.

The yellow precipitate may be titrated with permanganate by the Emmerton method involving solution in ammonium hydroxide, reduction of the acidified with  $H_2SO_4$  solution in a Jones reductor and oxidation of the reduced molybdenum compound by ordinary KMnO<sub>4</sub> solution. This is much used, but offers no decided advantages over the methods given above.

For the analysis of special steels, Johnson or Demorest should be consulted.

## SULFUR.

The sulfur assay generally involves the oxidation of S to SO<sub>3</sub>, and the subsequent precipitation of  $-SO_4$  by Ba++. As some SO<sub>2</sub> is evolved by the combustion of ordinary illuminating gas, it is advisable to use only alcohol or gasoline flames when determining sulfur. The crucible may, however, be placed in a large square of asbestos board which has been punched to receive it; and the regular Bunsen used.

It is advisable to pour the barium solution into the sulfate solution as less soluble barium salts will thus contaminate the precipitate. Even if this is carried out, BaSO<sub>4</sub> will include (enclose) and occlude (attract and stick to) BaCl<sub>2</sub>, alkali chlorids and  $(NH_4)_2SO_4$  appreciably if there is much BaSO<sub>4</sub> formed. Hence it is advisable to have as little precipitate as accuracy will permit.

(NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> will be lost on ignition, but occluded as well as included BaCl<sub>2</sub> will make the precipitate heavier. Occluded nitrates may cause high results. High iron concentrations especially if present as Fe<sup>III</sup> will cause high results unless considerable HCl is present, on the other hand BaSO<sub>4</sub> is slightly soluble in excess HCl. Too high temperature during ignition will give BaO.

The precipitate of BaSO<sub>4</sub> comes down so fine that digestion on the steam bath for an hour is advisable to prevent the precipitate from running through the filter. Or, the precipitate may be allowed to settle over night at room temperature.

Sample No. Analysis of

Procedure For Iron Ores, No. 4.—Mix I gram of ore with about 5 grams Na<sub>2</sub>CO<sub>8</sub> and 0.5-1 gram NaNO<sub>8</sub> in a 30-gram platinum crucible by stirring and adding the filter paper wipings to the crucible. Heat to quiet fusion, remove from heat and chill the outside of the crucible by dipping in water. The contents are then removed into a beaker containing about 75 cubic centimeters warm water and the crucible and cover washed out into and removed from the beaker. The melt is dissolved by heating; dilute HCl is added carefully until the effervescence ceases, a few drops of alcohol added to precipitate Mn, and the solution is filtered. To the filtrate is added 25 cubic centimeters of a 5 per cent. solution of BaCl<sub>2</sub>, or better pour the hot solution into 25 cubic centimeters of a 5 per cent. BaCl, solution diluted to 100 cubic centimeters. After digesting on a steam bath or allowing to settle over night, the supernatant liquid is decanted through the filter, and the BaSO4 is filtered off and weighed. It is advisable to dry the precipitate, separate it from the filter onto glazed paper and ignite the filter separately by winding around it a piece of platinum wire the end of which is sealed into a piece of glass tubing: holding the burning paper above the weighed crucible. The bulk of the precipitate is then brushed from the glazed paper into the crucible, and the latter ignited, cooled and weighed. It is also advisable to give the crucible an H<sub>2</sub>F<sub>2</sub> treatment (adding two drops of concentrated H<sub>2</sub>SO<sub>4</sub>) to rid of silica. (Since iron ores are apt to be very low in S, 2-5 grams of sample may be taken.)

The oxidation of sulfur in iron and steel is accomplished by not too rapid solution of 2-5 grams in 25-40 cubic centimeters concentrated HNO<sub>3</sub> in a platinum dish with subsequent addition of a few cubic centimeters of concentrated HCl and 0.5 gram KClO<sub>3</sub>. Evaporate to dryness, bake, dissolve in HCl and rebake. The Na<sub>2</sub>SO<sub>4</sub> may then be washed out of the Fe<sub>2</sub>O<sub>3</sub> by warm water, and the solution of Na<sub>2</sub>SO<sub>4</sub> poured into 5 per cent. BaCl<sub>2</sub> solution as above. Or the Fe<sub>2</sub>O<sub>3</sub> may be dissolved in HCl and the BaCl<sub>2</sub> solution added to this. In this latter case, it is advisable to reduce the iron to the ferrous condition by Al powder before precipitation of BaSO<sub>4</sub> as less contamination of the precipitate thus results.

Analysis of Sample No. Sulphur in 1.3736 grams pyrites ore may be oxidized by 10 cubic centimeters of a mixture of 2 parts by volume liquid Br and 3 parts CCl<sub>4</sub> in the cold. For mattes, concentrated HNO<sub>8</sub> saturated with KClO<sub>8</sub> may be used. The subsequent procedure is the same as in the case of iron and steel.

Sulfur may be determined volumetrically by evolution as H<sub>2</sub>S which is absorbed in ammoniacal CdBr<sub>2</sub> solution. This is rendered acid with HCl and the H<sub>2</sub>S changed to HI with iodine solution which has been standardized against bichromate solution of known strength. The end point is the blue color imparted to starch solution by free iodin.

#### CARBON.

The direct combustion method, using an electric furnace and an oxygen tank fitted with a diaphragm valve, is the method most often used for the determination of carbon in steel and iron. It is often superseded by metallographic methods.

Oxygen is purified by passing through KOH followed by H<sub>2</sub>SO<sub>4</sub>, Zn and CaCl<sub>2</sub> bottles and tubes in the order named. It then passes through a ½-inch silica combustion tube set in the furnace. The sample is placed in a groove made in a sand or alundum lining of a boat of platinum, silica, alundum or porcelain. The boat is put into the furnace through the end nearer the oxygen tank. It almost touches a larger boat containing the wire form of CuO which has been recently ignited. Several grams of PbCrO<sub>4</sub> should be added to this or a following boat. The ends of the combustion tube should project a few inches from the outer wall of the furnace and should be covered by wicks of blotting paper which dip into suspended beakers ¾ full of water.

The mixture of oxygen and carbon dioxid passes from the silica tube through 15-20 cubic centimeters of granulated Zn, and about the same amount of CaCl<sub>2</sub> into the weighing bulb. This contains 30 per cent. KOH solution in a series of bulbs which are manufactured in a number of forms. Those having sealed on CaCl<sub>2</sub> tubes are convenient, as some sort of weighable CaCl<sub>2</sub> tube must follow the bulb tube. The former is protected by a second CaCl<sub>2</sub> tube. Finally, an aspirator may be added to further regulate the flow of oxygen through the apparatus.

Procedure No. 5.—Weigh 2 grams of sample, place in a fairly compact mass in the groove of the boat, and after sweeping out the tube and weighing bulbs with oxygen for 15 minutes (the furnace being at 1,000° C.), place the boat in

Analysis of Sample No. the furnace as soon as the *weighed* bulbs have been replaced in position. In practice, the final weighing for one analysis serves as the initial weighing for the second, etc. Change in weight in KOH bulb  $\times$  0.2727  $\div$  2 = per cent. C in steel. It is hardly necessary to use rubber stoppers for the ends of the absorption vessels during weighing, especially if the latter have restricted openings.

In addition to the above, the color method is often used. 0.2 gram of two steels of similar manufacture and character, the carbon content of one of which is known, are dissolved in equal volumes of C. P. HNO<sub>2</sub> and the solutions are poured into two graduated colorimeter tubes. The unknown solution is then diluted until it is the same shade of brown or vice versa and the per cent. carbon calculated in terms of the water added. In practice the unknown is compared to a scale of known in comparison or test tubes.

### MANGANESE.

The bismuthate method is convenient for the determination of large or small quantities of manganese. It is based on the fact that HMnO<sub>4</sub> is stable in HNO<sub>8</sub> specific gravity 1.135 in the absence of -Cl. The method is generally applicable with the following precautions. -Cl must be rid of by boiling with HNO<sub>3</sub>. WO<sub>8</sub> of W steels must be filtered off after precipitation twice by evaporation to syrupy or dryness with aqua regia, and the latter replaced by HNO<sub>3</sub>. Va and Cr must be precipitated by CdCO<sub>3</sub> from a dilute H<sub>2</sub>SO<sub>4</sub> solution, and the latter replaced by HNO<sub>3</sub> by boiling with concentrated acid and diluting.

One gram of sample is usually taken; solution in dilute HNO<sub>8</sub> 1-3 is used for steels, H<sub>2</sub>F<sub>2</sub> is added to silicious materials, and ores are fused with 10 grams NaHSO<sub>4</sub>, 1 gram Na<sub>2</sub>SO<sub>8</sub>, and 0.5 gram Na<sub>2</sub>F<sub>2</sub>. An ether separation is advisable for alloy steels. (p. 24.)

Procedure No. 6.—Fuse or dissolve the sample, cool, add 0.5 gram of sodium bismuthate to the solution and heat till the pink color disappears. Clear the solution with sulfurous acid, heat to boiling to rid of excess of SO<sub>2</sub> and cool.

Add about 3 grams bismuthate and agitate for a few minutes, but do not heat. Add 50 cubic centimeters water (containing 30 cubic centimeters HNO<sub>3</sub> specific gravity 1.42 per liter), and filter off the bismuthate by running onto an asbestos or better a large alundum filter, and washing carefully. It is very convenient to

Sample No. Analysis of

use a one or two hole bell jar, and to suck the filtrate directly into the beaker. Run into the beaker a measured amount of standardized FeSO<sub>4</sub> solution and titrate the excess with KMnO<sub>4</sub> solution (1 gram per liter) which has been standardized best against a known Mn ore or steel.

FeSO<sub>4</sub> solution per liter: 12.4 grams FeSO<sub>4</sub>, 25 cubic centimeters H<sub>2</sub>SO<sub>4</sub> specific gravity 1.84, 25 cubic centimeters glacial phosphoric acid.

# CONDITIONS FOR THE TITRATION OF IRON BY PERMANGANATE.

Iron is usually titrated with KMnO<sub>4</sub>\* in a chlorid solution. (One gram of ore\*\* is dissolved in 15 cubic centimeters concentrated HCl, diluted, temperature being kept below 100° C. as FeCl<sub>3</sub> is volatile; and unless the ore contains Ti, the solution is run through a reductor.) This titration of Fe<sup>II</sup>—Fe<sup>III</sup> is subject to error from the following sources:

Improper standardization. Side reactions of KMnO<sub>4</sub>. Tardy action of KMnO<sub>4</sub>.

Iron wire is not thoroughly satisfactory for standardizing as the carbon in it forms hydrocarbons which reacting with the KMnO<sub>4</sub> give results varying from 99.8-100.6 per cent. Fe.

Instead of Fe wire either a known iron ore (a known steel or a known solution of Mohr salt) or best of all 0.25-0.3 gram of Na<sub>2</sub>C<sub>2</sub>O<sub>4</sub> of known quality may be used.

Procedure: dissolve some 0.3 gram of Na<sub>2</sub>C<sub>2</sub>O<sub>4</sub> in 200-250 of hot water (80-90° C.) add 10 cubic centimeters concentrated H<sub>2</sub>SO<sub>4</sub> 1-1, titrate with 0.560 per cent. KMnO<sub>4</sub> to first trace of pink. The first 0.5 cubic centimeter and the last 0.5-0.8 cubic centimeter of solution should be added dropwise, and each drop allowed to decolorize before adding the next. Stir vigorously and continuously.

Burette delivery should be from 10-15 cubic centimeters per minute. Final temperature should not be below 60° C.

<sup>\*</sup> Test problem.

<sup>\*\*</sup> A moisture test of the ore should be taken.

Analysis of	Sample No.		
		-	
	•		
<u> </u>		·	
			•
			<del></del>

The KMnO<sub>4</sub> solution should have been filtered through asbestos and kept away from the light, and free from organic matter. It should be standardized frequently.

A gram of NaOH per liter may be added to insure stability.

When standardizing against a known ore, the conditions are more nearly those of actual titration of (1 gram) unknown samples. The following conditions must be observed:

Temperature must be at least 70° C.

The concentration of HCl must be not more than 60 cubic centimeters dilute acid specific gravity 1.04 per 150 cubic centimeters of solution, i. e., 10 cubic centimeters HCl specific gravity 1.20 per titration.

The volume of the solution must be about 400 cubic centimeters per sample of about 0.4 gram Fe.

The beaker in which the titration is made should set on a white surface.

One gram of MnSO<sub>4</sub> must be present else side reactions due to HOCl formation by KMnO<sub>4</sub> will cause high results.

The addition of H<sub>2</sub>PO<sub>4</sub> (glacial, specific gravity 1.75) forms an almost colorless phosphate with the Fe and thus makes the end point more distinct.

Fifty cubic centimeters of the following solution is a convenient amount. 180 grams MnSO<sub>4</sub> dissolved in 1,300 cubic centimeters water; 350 cubic centimeters H<sub>2</sub>SO<sub>4</sub> specific gravity 1.84; 350 cubic centimeters H<sub>3</sub>PO<sub>4</sub> specific gravity 1.75. By observing these precautions, HOCl will not be formed if the titration is carried on rapidly in warm acid solution. In cold or non-acid solution the reaction will not run quantitatively; *i. e.*, instead of

$${\rm 10FeSO_4 + 2KMnO_4 + 8H_2SO_4 = 5Fe_1(SO_4)_3 + K_2SO_4 + 2MuSO_4 + 8H_2O_4 + 8H_$$

we find that the MnO<sub>2</sub> of the KMnO<sub>4</sub> will either oxidize some of the HCl or it will precipitate as MnO<sub>2</sub>.

Upon solution, all the iron is in the ferric state. It is reduced to ferrous iron by adding 20 mesh Zn to the solution, and allowing it to dissolve completely, or better still, by sucking the solution through a column of 20 mesh Zn (previously

Analysis of Sample No. amalgamated by stirring in a beaker with a 1 per cent.  $HgCl_2$  solution containing a globule of mercury) weighing about  $\frac{2}{3}$ -pound contained in a 1-inch glass tube. This tube (the Jones Reductor) should have a stop-cock at the base and be kept full of wash solution, when not in use. Wash solution = 50 cubic centimeters  $H_2SO_4$  specific gravity 1.84 per 2 liters water.

Solutions of FeCl<sub>a</sub> should not be boiled as FeCl<sub>a</sub> is appreciably volatile with steam.

If the end point is run by, add I cubic centimeter FeSO<sub>4</sub> solution, titrate to pink, read, add another cubic centimeter FeSO<sub>4</sub>, titrate to pink again, read, and deduct the difference between the two readings from the first reading.

Ores containing appreciable amounts of Ti cannot be reduced with Zn as the latter reduces the Ti also. The Ti is then oxidized by the KMnO<sub>4</sub> causing high results. Such ores should be reduced by SnCl<sub>2</sub>.

o.88 per cent.  $K_2Cr_2O_7$  solution may also be used for determining ferrous iron.  $K_2Cr_2O_7$  is used very little compared to KMnO<sub>4</sub>. In this procedure the iron is dissolved in HCl, SnCl<sub>2</sub> solution is added during solution until finally a colorless (cooled) solution is obtained. The slight excess of SnCl<sub>2</sub> in the cooled solution is oxidized with HgCl<sub>2</sub> solution, and the iron titrated with either KMnO<sub>4</sub> or more usually with  $K_2Cr_2O_7$ . This latter procedure involves the use of a tile covered with drops of very dilute fresh  $K_3Fe(CN)_3$  solution. When the end point is approached drops of the titrated solution touched to the indicator with the tip of a stirring rod will give a faint and finally no blue color. This is the end point. Diphenylcarbazid may be used as an inside indicator. A potentiometer gives an excellent end point. (This may be done as Optional Procedure No. A.)\*

The insoluble residue from the HCl solution should always be filtered off; washed with hot dilute HCl and fused with Na<sub>4</sub>CO<sub>5</sub> or KHSO<sub>4</sub>. Even colorless residues may contain Fe as phosphate. The Fe thus dissolved is added to the main solution. If none is found the residue is rejected.

Care must be taken in washing as cold or non-acid water allows the formation of insoluble Fe salts in the filter paper.

Analysis of	Sample No.		•	
· · · · · · · · · · · · · · · · · · ·	······································			
			-	
_				
	•			
			•	
· -				
	<del></del>		<del></del>	
		,		
,				
- v. u				
				]

.

# IRON ORES, IRON AND STEEL, LESSER CONSTITUENTS.

## NICKEL.

The dimethylglyoxime method for nickel is very accurate, and when the precipitate is dissolved and titrated, the method has the added advantage of rapidity. Furthermore, it is applicable to practically all steels containing nickel, to nichrome, to monel metal, to brasses and to German silver.

Procedure No. 7.—Dissolve in a large beaker I gram drillings in 20 cubic centimeters HCl I-I and add IO cubic centimeters dilute HNO<sub>3</sub>. Boil off excess of acids and add 30 cubic centimeters H<sub>2</sub>SO<sub>4</sub> containing a little HNO<sub>3</sub>. Add IO-I5 cubic centimeters powdered tartaric acid (citric acid may also be used), when all dissolved, make slightly alkaline with NH<sub>4</sub>OH. Then add HAc until the solution reacts acid (if Ni only is wanted the most of the Fe may well be separated by ether from HCl solution specific gravity I.IOO-I.IO5). Where-upon heat to boiling and add 5 times as much dimethylglyoxime by weight as I per cent. solution in alcohol as there is Ni in the sample. Make the solution alkaline with NH<sub>4</sub>OH, and if convenient allow the precipitate formed to stand an hour before filtering into a gooch. There should be 4-5 times as much solution as there is precipitating solution added to prevent the solution of the latter and the precipitate in the former. The Ni precipitate should not be allowed to get cold.

In routine work, the precipitate may be filtered off at once into a weighed alundum or gooch crucible, washed with very dilute NH<sub>4</sub>OH and with a little alcohol and weighed. Precipitation is 20.32 per cent. Ni.

The analysis may be further hastened by catching the precipitate on alundum or on a Büchner funnel covered with an asbestos mat, dissolving the precipitate in a little warm concentrated HNO<sub>8</sub> catching the solution (using a bell jar) in a beaker, boiling with a few grams (NH<sub>4</sub>)<sub>2</sub>S<sub>2</sub>O<sub>8</sub> and cooling with clean pieces of ice. The indicator is 10 cubic centimeters of 0.05 per cent. AgNO<sub>3</sub> plus 10 cubic centimeters 2 per cent. KI solution, and the titrating solution, 0.475 per cent. KCN solution which has been standardized against a steel, the Ni content of which is known.

Analysis of Sample No.

# COPPER."

Procedure No. 8.—3-5 grams of steel are dissolved in 7 cubic centimeters  $H_2SO_4$  diluted with 60 cubic centimeters water. Drop into the solution a piece of sheet Al 1½-inch by 1½-inch with turned over edges so that both sides will be exposed. Boil the solution 25 minutes to precipitate the Cu, and decant the solution through a filter, washing the Cu in the beaker with pure water. Rinse off the Al with dilute HNO<sub>3</sub> and remove with a glass rod; place the beaker under the washed filter, poke a hole in the latter, and wash any adhering Cu into the beaker. Electrolyze the dissolved Cu, depositing on a gauze cathode previously weighed. 0.5-gram Cu will deposit quantitatively from a nitrate solution containing a few drops only excess of HNO<sub>3</sub> in from 15 minutes to an hour if enough surface is provided. Add a few cubic centimeters  $H_2O_2$  solution towards the end of the deposition to rid of HNO<sub>2</sub> in which Cu is somewhat soluble.

A 6-volt 100-ampere battery is a convenient source of current, a larger one is better; 110 volts direct current may also be used if suitable resistance is inserted in the circuit. For rapid determinations, using more than a small fraction of an ampere, plain electrodes and motionless solutions must be superceded by one of the following: Rotating the anode or cathode, rotating the solution by means of a small air jet or a solenoid, or increasing the surface of the electrodes by using Pt. (or Cu plated brass) gauze. If the gauze be used, there is not much saving in time in using rotating devices. A rheostat and meters are advisable.

The procedure details are to weigh the clean electrodes after washing it in distilled water followed by alcohol (and ether) and air blast or air bath drying; insert in the solution, electrolyze, decrease current, withdraw the electrode, while rinsing it with water, set in a beaker of pure water a minute, withdraw, rinse with alcohol (and ether), dry and weigh. Difference in weight is Cu.

Instead of precipitating with Al, H<sub>2</sub>S may be used. The precipitate is filtered through a small paper washed with warm dilute H<sub>2</sub>SO<sub>4</sub>, ignited in porcelain, dissolved in HNO<sub>3</sub>, transferred to a (weighed) Pt crucible, treated with H<sub>2</sub>F<sub>2</sub> and H<sub>2</sub>SO<sub>4</sub> and reweighed or dissolved in H<sub>2</sub>SO<sub>4</sub> or HKSO<sub>4</sub> (warmed) and electrolyzed or determined colorimetrically. If Mo is present this last solution must be treated with NaOH to separate the copper as CuO. The CuO is filtered off and the Cu determined as above. The CuO is filtered from the filtrate as sulfid and roasted to oxid.

### TUNGSTEN."

**Procedure No. 9.**—Two grams of the steel\* may be dissolved in dilute HCl or H<sub>2</sub>SO<sub>4</sub> in a casserole, and the WO<sub>3</sub> precipitated by the addition of 60 cubic

<sup>\*</sup>W steels may be softened by annealing at 700° for 2 hours.

Analysis of	Sample No.		
		 <del> </del>	
•			
,			
		<del></del>	

.

centimeters HNO<sub>3</sub> 1.42. Digest the WO<sub>3</sub> until it is of a clear yellow color: to the solution a little cinchonine solution\* is added if W only is wanted, (or else the solution is taken down to dryness. Redissolve in HCl 1-1 to which 10 cubic centimeters concentrated HNO<sub>3</sub> is added), dilute and filter. (Dissolve and adhering WO<sub>3</sub> to casserole with a bit of filter paper dipped in NH<sub>4</sub>OH and add to main precipitate) wash precipitate with water place in weighed platinum crucible and ignite, cool, weigh, treat with H<sub>2</sub>F<sub>3</sub> and H<sub>2</sub>SO<sub>4</sub>, reignite on a ring burner, cool and weigh. WO<sub>3</sub> is volatile at bright red heat, hence both the above ignitions should be done carefully WO<sub>3</sub>  $\times$  79.31  $\div$  2 = W. (Fusion with HKSO<sub>4</sub> or Na<sub>2</sub>CO<sub>3</sub> is necessary to rid the precipitate quantitatively of iron.)

The WO<sub>3</sub> being an acid anhydrid may be dissolved in NaOH and the excess of the latter titrated with standard HCl using phenolphthalein as an indicator, and standardizing against a steel the W content of which is known.

## CHROMIUM ONLY."

(Chromium may be determined by treating the HCl solution of 2 grams of steel as follows: The solution is evaporated to specific gravity 1.100-1.105 and the iron separated (twice), by shaking in a separatory funnel, with ether. Add NH<sub>4</sub>OH to the water layer, filter, dry ignite and fuse with Na<sub>2</sub>CO<sub>3</sub> and NaNO<sub>3</sub>. Extract with water, and wash with dilute Na<sub>2</sub>CO<sub>3</sub> solution giving Na<sub>2</sub>CrO<sub>4</sub> and some Al in the filtrate. Evaporate to syrupy with excess NH<sub>4</sub>NO<sub>3</sub>. Dilute and filter again, giving Na<sub>2</sub>CrO<sub>4</sub> only in the filtrate. Boil to rid of NH<sub>4</sub>OH and reduce the Cr with H<sub>2</sub>SO<sub>3</sub>. Boil and add NH<sub>4</sub>OH in excess, filter, wash, ignite and weigh the Cr<sub>2</sub>O<sub>3</sub> 68.46 per cent. Cr.)

## VANADIUM ONLY."

For Ores: Dissolve I gram of ore in 25 cubic centimeters of the following: 300 cubic centimeters HNO<sub>8</sub>, specific gravity 1.42; 300 cubic centimeters H<sub>2</sub>SO<sub>4</sub>, 1.84; 300 cubic centimeters water 100 cubic centimeters H<sub>3</sub>PO<sub>4</sub> specific gravity 1.75, 4-5 grams MnSO<sub>4</sub>; dilute and filter. To filtrate, add 1-gram sodium bismuthate, stir, heat to rapid boiling. Add a few cubic centimeters HCl and re-boil a moment. Cool, ice may be used, filter. Add a pipetteful of standardized

<sup>\*</sup>Twenty-five grams cinchonine in 200 cubic centimeters dilute HCl 1-1.

Sample No. Analysis of

ferrous sulfate solution, and titrate with KMnO<sub>4</sub>, running a known sample at the same time.

For Steels: No. 161 Bureau of Standards method involves the separation of Fe by absorption in a Hg cathode. (See p. 30.) The Cr is absorbed by the Hg, hence it cannot be determined in the same sample. Cr should be determined by Cain's method<sup>12</sup> as PbCrO<sub>4</sub>, although the combined method, given below, if practiced becomes accurate enough for routine work.

## CHROMIUM AND VANADIUM IN STEEL. USUAL METHOD."

Procedure No. 10.—Dissolve 2 grams of known and unknown steel (or 0.5-gram ferro) in 10 cubic centimeters H<sub>2</sub>SO<sub>4</sub> and 50 cubic centimeters water in a large beaker. Oxidize the solution by adding 60 cubic centimeters dilute HNO<sub>2</sub> '1-1 and boil to rid of nitrous fumes. Dilute to about 300 cubic centimeters and add enough KMnO<sub>4</sub> solution (30 grams KMnO<sub>4</sub> per liter) to keep precipitated a little MnO<sub>2</sub> after boiling 20 minutes. Filter off the MnO<sub>2</sub> and wash with very dilute H<sub>2</sub>SO<sub>4</sub>.

To filtrate, add standardized n/50 FeSO<sub>4</sub> solution until the color of the solution changes to colorless or light green, add a few tenths of a cubic centimeter more, and titrate the excess with n/50 KMnO<sub>4</sub>. (Fe equivalent of FeSO<sub>4</sub> — Fe equivalent of KMnO<sub>4</sub>)  $\times$  0.3103 = amount Cr. Better compare direct to standard.

Va. Add exactly 1 cubic centimeter of a 0.1 per cent. solution of  $K_s$ Fe  $(CN)_6$   $_3H_2O$  to the beaker and add the standard  $FeSO_4$  until the yellow shade of the Cr changes to a full green. Do not titrate to blue, but titrate the known and unknown to the same color or Fe equivalent  $FeSO_4 \times 0.914 =$  amount Va. Better compare to standard. These colors are difficult to distinguish without practice.

#### **MOLYBDENUM.**"

**Procedure No. 11.**—Dissolve 2 grams of steel in 50 cubic centimeters HCl, heat to boiling and add a little HNO<sub>3</sub> to complete the solution and oxidize the Fe,

Sample No. Analysis of

but do not add excess. Evaporate to syrupy, add 50 cubic centimeters dilute HCl 1-3 filter off WO<sub>2</sub> if present washing with dilute HCl. If much WO<sub>2</sub> is present a second evaporation is advisable.

To the filtrate add dilute NaOH to neutralize nearly all the HCl, but do not add enough to cause the color of the solution to deepen. Transfer to a separatory funnel, and drop the solution into a 500 cubic centimeter graduated flask containing 150 cubic centimeters hot NaOH solution, 6 per cent., shaking the flask vigorously while dropping in the NaOH so that the Fe(OH)<sub>3</sub> formed will not carry down with it any Mo. Wash out the separatory funnel into the flask.

Dilute the solution to the 500 cubic centimeter mark, and filter through a large paper washing well. To 250 cubic centimeters of the filtrate, add HCl to acidity (methyl orange), add a few cubic centimeters more, followed by a few drops sulfurous acid to reduce Cr, and boil.

Add 10 grams AmAc stir to dissolve, add 40 cubic centimeters of a 1 per cent.  $PbAc_2$  solution, filter, dry the precipitate, separate, ignite paper, and bulk of precipitate and ignite both to faint redness in a weighed porcelain crucible. Two  $\times$  26.16 per cent. of precipitate is Mo.

A somewhat shorter method is to precipitate the Mo as MoS from dilute H<sub>2</sub>SO<sub>4</sub> solution, filter through asbestos, wash, dissolve the precipitate in aqua regia, and evaporate the solution with H<sub>2</sub>SO<sub>4</sub> (if Fe remains the H<sub>2</sub>S treatment should be repeated), rid of Cu, if present with NaOH precipitate the Mo as MoS filter, wash and cautiously ignite to MO<sub>3</sub>.

## TITANIUM."

Procedure No. 12.—Dissolve 5 grams of known and unknown sample of iron in 40 cubic centimeters HCl specific gravity 1.05, dilute, filter and wash. The insoluble contains all the Ti. Ignite the paper and precipitate in a 20-gram Pt crucible and treat with H<sub>2</sub>F<sub>2</sub> and H<sub>2</sub>SO<sub>4</sub>. Fuse the residue with Na<sub>2</sub>CO<sub>3</sub> dissolve the fusion in water, filter and dissolve the insoluble in H<sub>2</sub>SO<sub>4</sub> (if necessary HKSO<sub>4</sub>) and combine the two solutions.

Add H<sub>2</sub>O<sub>2</sub> to both the known and the unknown solutions and compare in a

Analysis of Sample No. colorimeter. If a Fe free Ti solution is used for a standard, ferric alum should be added to it to match the iron color of the unknown, before H<sub>2</sub>O<sub>2</sub> is added to either.

Ores may be fused with Na<sub>2</sub>CO<sub>3</sub> the washed residue after dilute HCl treatment fused with HKSO<sub>4</sub> and Ti determined colorimetrically in the latter dissolved fusion. Ferros are best determined by one of the gravimetric processes.

## **ALUMINIUM.**"

Procedure No. 13.—The chlorid solution of 5 grams of steel is evaporated with H<sub>2</sub>SO<sub>4</sub> to small bulk, diluted to 200 cubic centimeters and almost neutralized with NH<sub>4</sub>OH. The solution is then poured into a medium sized beaker containing a layer of about 300 grams Hg and a contact in the form of a sealed in Pt wire at the base, (or else a glass tube passing down through the solution, with a sealed in Pt wire at the lower end, containing either Hg of a place of Cu wire attached to the inside end of the Pt wire), forming the electrical contact between the current supply and the Hg cathode.

The Fe is driven in to the cathode (with a rotating anode a current of several amperes may be used and the Fe rid of in an hour) leaving Al, P, and Ti only in solution. Precipitate twice with NH<sub>4</sub>OH and ignite and weigh in a porcelain crucible. The P<sub>2</sub>O<sub>5</sub> should be determined in another sample, the above residue examined colorimetrically for Ti as TiO<sub>2</sub>, and these deducted giving total Al<sub>2</sub>O<sub>3</sub>.

Peters method is not recommended for the determination of Al<sub>2</sub>O<sub>3</sub>. The latter may be obtained by a phenylhydrazone precipitation, by an ether separation followed by acetate separation and ammonia precipitation, or finally the electrolytic procedure may be used.

Analysis of	Sample No.	
A 1		
		-

.

1

## LIMESTONES, SLAGS AND CLAYS."

## LIMESTONES.

Procedure No. 14—Insoluble Silicious Matter.—Boil I gram of limestone in a beaker with 5 cubic centimeters HCl and I cubic centimeter HNO<sub>3</sub> plus 25 cubic centimeters H<sub>2</sub>O until action ceases. Wash off cover glass into beaker and evaporate solution to dryness. Redissolve in 10 cubic centimeters concentrated HCl diluted with 50 cubic centimeters H<sub>2</sub>O. Filter. Save filtrate; wash precipitate well with hot water, ignite and weigh. Examine for metals. If SiO<sub>2</sub> is wanted, baking, followed by solution and H<sub>2</sub>F<sub>2</sub> treatment are necessary.

Al<sub>2</sub>O<sub>8</sub> and Fe<sub>2</sub>O<sub>8</sub>. Heat the filtrate to boiling, add slight excess NH<sub>4</sub>OH. Boil a few minutes, filter and wash. Dissolve precipitate in a little HCl on the paper, wash and let the always slightly acid solution run into a small beaker. Reprecipitate with slight excess NH<sub>4</sub>OH, filter, wash, ignite and weigh as Fe<sub>2</sub>O<sub>3</sub>; Al<sub>2</sub>O<sub>3</sub>. Fe may be titrated with KMnO<sub>4</sub> after running the HCl solution through a reductor.

Ca0. Unite the two filtrates, heat to boiling, add 25 cubic centimeters 5 per cent.  $(NH_4)_2$   $C_2O_4$  solution. Digest on the steam bath as long as convenient up to 1 hour. Filter, wash with hot water. Either ignite to CaO in a weighed Pt crucible or else after washing the precipitate into 100 cubic centimeters warm water, by poking a hole in the filter paper; and adding 30 cubic centimeters N/10  $H_2SO_4$ , to the filtrate to form  $H_2C_2O_4$ , titrate with KMnO<sub>4</sub>.

Mg0. To the above filtrate add 30 cubic centimeters saturated solution microcosmic salt. Acidulate with HCl, and evaporate to not more than 300 cubic centimeters. Cool the still acid solution. Add NH<sub>4</sub>OH to alkalinity (litmus). Add 5 cubic centimeters NH<sub>4</sub>OH, and digest on the steam bath, or stand over night. Filter and wash with dilute NH<sub>4</sub>OH. Dry and separate from paper onto glazed paper, burn filter paper on a platinum wire over a weighed porcelain crucible; combine the precipitates and weigh as Mg<sub>2</sub>P<sub>2</sub>O<sub>7</sub>. 36.22 per cent. MgO.

A moisture test at 100° C. and a loss on ignition test at the temperature of the blast lamp should be made on separate 1 gram specimens of the limestone.

Analysis of	Sample No.	
		·
	· · · · · · · · · · · · · · · · · · ·	
	,	

•

Ca can generally be determined by solution in HCl, evaporation to dryness (to separate  $SiO_2$ ) re-solution in HCl, precipitation of Fe and Al by  $NH_4OH$ ; and precipitation of  $Ca^{++}$  with  $--C_2O_4$ . This may be weighed or titrated.

Ca and Mg may be determined in the same solution as follows:<sup>18</sup> Precipitate the SiO<sub>2</sub>, Fe and Al as above using a slight excess only of NH<sub>4</sub>OH. To the filtrate, add excess H<sub>2</sub>C<sub>2</sub>O<sub>4</sub> solution followed by dilute NH<sub>4</sub>OH to alkalinity (litmus). Add enough ammonium arsenate to precipitate the Mg and add NH<sub>4</sub>OH constantly stirring the hot solution until the magnesium ammonium arsenate comes down; *i. e.*, about 10 cubic centimeters concentrated NH<sub>4</sub>OH. Cool add 10 cubic centimeters more NH<sub>4</sub>OH.

Ca is determined as follows: The ammonium arsenate solution should be added right after the oxalic acid so that Mg may not precipitate as oxalate. Let stand 20 minutes (up to 12 hours) and filter on paper and wash with dilute ammonia water. Transfer washed precipitate to ½ or ½-liter Erlenmeyer, washing with hot water followed by 1 per cent. H<sub>2</sub>SO<sub>4</sub> to dissolve. Make up the volume in the flask to 74.80 cubic centimeters and add 10 cubic centimeters 1-1 H<sub>2</sub>SO<sub>4</sub>. Heat to 75° C., and titrate with 0.566 per cent. KMnO<sub>4</sub> solution. Standardized with Ca and Mg precipitated under the same conditions.

Mg. Cool the Erlenmeyer in which the titration was made, add 25 cubic centimeters I-I H<sub>2</sub>SO<sub>4</sub>, bulk equalling about 125 cubic centimeters; add 5 grams powdered Kl slowly and immediately titrate to straw color with 1.9 per cent. Na<sub>2</sub>S<sub>2</sub>O<sub>2</sub> solution, titration taking 8 minutes. Standardize the Na<sub>2</sub>SO<sub>3</sub> against MgO, Mg ribbon, or best against a known limestone as above.

Finally, the Mg may be determined volumetrically using a milli-voltmeter as an indicator. This procedure is extremely rapid and valuable. It should be studied by the reader, but attempted only by those familiar with elementary electrochemistry.

Optional Procedure No. B.—The HCl solution of the limestone<sup>10</sup> is titrated with NaOH and the constancy of concentration of [H]+ or -[OH], (measured on the millivoltmeter and plotted as ordinate) due to the precipitation of Mg(OH)<sub>2</sub> is plotted on paper or automatically, against the standardized alkali added (as abissa).

## SLAGS.

Si0<sub>2</sub>. See p. 3.

Cu. Disintegrate with HCl followed by HNO<sub>3</sub>. Dehydrate SiO<sub>2</sub> with H<sub>2</sub>SO<sub>4</sub> boiling solution to fumes, cool, dilute and filter off. Neutralize H<sub>2</sub>SO<sub>4</sub> with NH<sub>4</sub>OH to only slight acidity, and throw down the Cu with Al. This latter should be Fe free or else the Fe in it determined carefully.

Fe. The Fe if present is, after reduction in the filtrate as Fe<sup>II</sup>. This may be titrated or oxidized with Br or HNO<sub>3</sub> and precipitated with NH<sub>4</sub>OH.

Analysis of Sample No. Ca. Ca may be determined in the last filtrate, but is more usually determined as follows. The solution freed from SiO<sub>2</sub> and oxidized with Br or HNO<sub>3</sub> add NH<sub>4</sub>OH to the first trace of Fe(OH)<sub>2</sub>. Add only enough hot saturated H<sub>2</sub>C<sub>2</sub>O<sub>4</sub> solution to dissolve this, and add NH<sub>4</sub>OH in slight excess. Add dilute warm H<sub>2</sub>C<sub>2</sub>O<sub>4</sub> to just dissolve the Fe(OH)<sub>2</sub> leaving the precipitate of CaC<sub>2</sub>O<sub>4</sub> in a clear yellowish green solution. Digest as long as convenient, decant through the filter, wash precipitate, and dissolve into a beaker with warm dilute H<sub>2</sub>SO<sub>4</sub>. Dilute to 200 cubic centimeters add 30 cubic centimeters, 10 per cent. H<sub>2</sub>SO<sub>4</sub>, and titrate with 0.569 per cent. KMnO<sub>4</sub> solution standardized against a known slag.

## Na2O AND K2O IN CLAYS.

Procedure No. 15.—The J. Lawrence Smith method for clays. Mix 1-gram clay with 1-gram NH<sub>4</sub>Cl in a porcelain mortar, grinding thoroughly. Add 8 grams of powdered CaCO<sub>3</sub> alkali free, and regrind. Pour into a 40 gram Pt crucible, lined with a little powdered CaCO<sub>3</sub> and sweep out the mortar with a little more CaCO<sub>3</sub> adding this to the charge. Place the crucible  $\frac{3}{4}$  through a piece of punched asbestos board, and put on cover.

Heat the crucible gently with the tip of a Bunsen flame until the NH<sub>2</sub> is all driven off. Then heat for 34-hour at full red heat. Cool, the melt and transfer quantitatively, washing with water, to a casserole. Digest the charge thoroughly, filter and wash with hot water until the filtrate amounts to 250 cubic centimeters.

To the hot filtrate, containing all the alkalis as chlorids, add NH<sub>4</sub>OH and (NH<sub>4</sub>)CO<sub>3</sub> in excess; digest a short while, and filter off the CaCO<sub>3</sub>. Wash latter with water containing a little NH<sub>4</sub>OH, concentrate the filtrate to small volume in silica or platinum, transfer to a Pt crucible, and ignite to dryness, using a ring burner. Weigh as NaCl+KCl. The chlorids should be white, and should leave no residue on solution in 10 cubic centimeters water. Add excess H<sub>2</sub>PtCl<sub>6</sub> solution, evaporate again to syrupy carefully. Add 20 cubic centimeters, 80 per cent. C<sub>2</sub>H<sub>5</sub>OH, and let stand until the Pt salts dissolve. Filter through a weighed gooch and wash the K<sub>2</sub>PtCl<sub>6</sub> with 80 per cent. C<sub>2</sub>H<sub>5</sub>OH, dry and weigh. K<sub>2</sub>PtCl<sub>6</sub> is 19.38 per cent. K<sub>2</sub>O. Eighty per cent. C<sub>2</sub>H<sub>5</sub>OH should be used and the Na<sub>2</sub>PtCl<sub>6</sub> given time to dissolve.

Analysis of	Sample No.	
	•	
	-	
	-	
1		

## NON FERROUS ALLOYS."

# BRASS, BRONZE, GERMAN SILVER, SOLDER, BEARING METALS.

Brass consists of Cu, Zn, Pb, Sn, (Fe, Al, Mn, P, Ni in small amounts).

Bronze consists of Cu, Sn, Zn, (P, Mn, Pb sometimes).

German Silver consists of Ni, Cu, Zn, Fe, Mn, (Co).

Solder consists of about half and half Pb and Sn, carrying a fraction of a per cent. of Sb, Zn, Cu and Bi.

Bearing metals consist of Pb, Sb, Cu, Sn, (Zn).

The above compositions are only approximate, and often overlap. Unusual constituents are often found, being added to impart certain specifications, or finding their way into the charge by accident.

#### BRASS.

Procedure No. 16.—Dissolve I gram of sample in 20 cubic centimeters HNO<sub>3</sub> I-I, and digest at about 90° C. for from 15-90 minutes. If the Sn content is less than 10 per cent. the hydrated SnO<sub>2</sub> formed\* will generally be as pure\*\* as the limit of accuracy of the entire analysis requires. Filter through paper or weighed alundum, wash with hot water containing a little HNO<sub>3</sub>, ignite and weigh as SnO<sub>2</sub> 78.81 per cent. Sn.

**Pb. Cu.** Electrolyze the filtrate, (if it contains no Mn,) after boiling off almost all the HNO<sub>3</sub>, and diluting to about 150-200 cubic centimeters. A Pt gauze anode and cathode are advisable. Rotation of the solution is convenient, but not necessary.\*\*\* As mentioned on p. 22 a difference of potential of about 6 volts and

<sup>\*</sup> Two to three cubic centimeters of a solution of egg albumin (white of one egg beaten with about 100 cubic centimeters water) will cause fine SnO<sub>2</sub> precipitates to coagulate.

<sup>\*\*</sup> The purification is described under bronze.

<sup>\*\*\*</sup> Current density is the ratio of the current-strength to the area of the electrode surfaces. A current-strength of one ampere per 100 square centimeters is the unit. The usual gauze cathode will allow  $I - 2 \times I$  the unit current density called N.D.100, but by rapid rotation up to 15 amperes may be used per 100 square centimeters of gauze surface, i. e., N.D.100 = 15. (Texts of Smith and Classen.)

Decomposition potential is the lowest potential which will allow a permanent plate of metal to deposit in a given solution. It is about 2 volts for most ordinary salts. Since such a low potential will not always allow a very large current to flow through the bath, it is customary to employ a battery of at least three lead storage cells in series.

Analysis of	Sample No.		
	-		
	·		
· · · · · · · · · · · · · · · · · · ·			
•			
		•	
			•

•

a current of 0.5-6 amperes is necessary. Add 1-2 cubic centimeters H<sub>2</sub>O<sub>2</sub> solution towards the end of the deposition.

The Pb deposits as PbO<sub>2</sub> 86.4 per cent. Pb on the anode. The anode is dipped in pure water, dried at about 230° C. and weighed. (PbO<sub>2</sub> is dissolved in HNO<sub>3</sub> to which H<sub>2</sub>O<sub>2</sub> solution, C<sub>2</sub>H<sub>5</sub>OH or H<sub>2</sub>C<sub>2</sub>O<sub>4</sub> solution has been added.) The Cu precipitates as metals on the cathode; it is washed with pure water dried with C<sub>2</sub>H<sub>5</sub>OH and (C<sub>2</sub>H<sub>5</sub>)O at 90° C. cooled and weighed. The electrolysis should take from 15-90 minutes depending on the current, electrodes and sample. Currents as high as 8-15 amperes have been used, and 0.5-gram Cu deposited in 10 minutes. Thirty-forty minutes is a conservative time for the deposition of the Cu in 1-gram brass using gauze electrodes, 2-4 amperes, rotation by means of a solenoid, motor or air blast.

Pb is also determined as sulphate, by boiling off the HNO<sub>3</sub> (in which PbSO<sub>4</sub> is slightly soluble), diluting the solution and pouring into dilute H<sub>2</sub>SO<sub>4</sub>. The same conditions as those for the BaSO<sub>4</sub> precipitation should be noted. Twenty cubic centimeters alcohol should be added to bring down the last of the PbSO<sub>4</sub>. PbSO<sub>4</sub> should not be ignited in platinum.

The following was given by the Bureau of Standards.

". . . the tin . . . (about I per cent. of brass sample No. 37) . . . contained impurities to the amount of 0.05-0.08 of I per cent., principally Fe, with some Cu, etc." (This is beyond the limit of accuracy of the Cu determination. Therefore the SnO<sub>2</sub> needs no further purification unless present in large amounts. Bureau uses I-I HNO<sub>3</sub>.)

"It is believed that the PbSO<sub>4</sub> gives the nearer actual Pb content. . . . the impurities in the PbO<sub>2</sub> are greater than the difference between the Pb obtained as PbSO<sub>4</sub> and the Pb obtained as PbO<sub>2</sub>."

When the Cu is all deposited from the solution (best determined by a calomel electrode and a mille-voltmeter),\* the electrodes are washed off into the solution and the following procedures gone through.

Fe. The sample is now examined for Fe by adding slight excess NH<sub>4</sub>OH filtering, and determining the washed precipitate as Fe<sub>2</sub>O<sub>8</sub> or dissolving, reducing and titrating with KMnO<sub>4</sub>.

<sup>\*</sup>Wait until you are sure all the Cu is deposited, and then wait three minutes more.

Analysis of Sample No. · Ni. Ni is determined as dimethylglyoxime. It is advisable to dissolve this and reprecipitate to clear of Zn. The conditions for the precipitation are given on p. 20. Tartaric acid need not be used. (Ni may be plated from a slightly alkaline with NH<sub>4</sub>OH solution, although Fe, if present, will also be carried down, by the current, and hence must first be separated.)

Zn. From the solution freed from Ni and Fe as above the Zn is separated by one of three methods.<sup>21</sup> Titration with K<sub>4</sub>Fe(CN)<sub>6</sub> solution, precipitation as phosphate, or deposition by the electric current. (Zn is sometimes precipitated as ZnS and roasted to oxid, but this compared with the above is rare.)

All determinations of Zn are best carried on in pure solution, subject to the following conditions.

Electrolysis. Boil off NH<sub>4</sub>OH with HNO<sub>3</sub> and rid of latter with H<sub>2</sub>SO<sub>4</sub> boiling to first fumes. Dilute to about 150-200 cubic centimeters and add 5-12 grams solid KOH to the solution. Electrolyze the heated complete solution using a solenoid or a rotating cathode and up to 5 amperes. A cathode of either Ni or Cu plated No. 100 brass gauze should be used. Platinum is not altogether suitable. The following precautions should be observed in Zn electrolysis.

All acids except H<sub>2</sub>SO<sub>4</sub> should be absent. Only a little excess H<sub>2</sub>SO<sub>4</sub> should be present, as alkali sulfates are not very soluble in alkaline solution.

Zn is lost during the washing off of the cathode with water. The formation of Zn(OH)<sub>2</sub> about makes up for this, if the washing is done rapidly and carefully. The procedure is to remove the electrode from the bath, dip into distilled water, and twirl rapidly between the fingers.

Drying is accomplished by dipping in  $C_2H_5OH$  then in  $(C_2H_5)_2O$  and drying above room temperature at about 80-90° C. over a hot plate. Else  $Zn(OH)_2$  will form.

The anode is best placed above the cathode. The latter or the solution should be rotated. Anode rotation is not satisfactory.

Time of electrolysis, about 40 minutes.

At end of 30 minutes, add ice to cool the solution. (As above, the addition

Analysis of Sample No. of the solid KOH heats the solution to about the right temperature for most of the electrolysis.)

Nitrates, chlorids and NH<sub>4</sub> salts should be absent.

The following remarks of Sand<sup>2a</sup> especially those on cooling, apply to zinc deposition by electrolysis. "The formation of crystals is almost always bound up with the intermediate formation of a metastable state. Thus crystallization from a solution is probably always connected either with local supersaturation or with the intermediate production of an undercooled liquid. In a similar way, it seems very probable that metal deposits are produced electrolytically in the first instance in a metastable state, and that it depends on the ease and rapidity with which this state is destroyed whether a fine or a coarse grained deposit will result."

Precipitation. The gravimetric procedure for Zn in pure solution is as follows: Add to the Zn solution as chlorid or nitrate, bulk about 150 cubic centimeters, NH<sub>4</sub>OH to slight alkalinity. Add a drop of litmus solution and drop in HNO<sub>4</sub> 1-5 till the color just changes to red. Exact neutrality must be reached.

Heat the solution and pour in  $(NH_4)_2HPO_4 = 12 \times$  the weight of Zn present. This is prepared from the commercial salt by adding to its solution,  $NH_4OH$  until a drop of phenolphthalein just turns pink. The object is to get exact neutrality, as here also the precipitate is soluble in both acid and alkali.

Digest the resulting precipitate on the steam bath for 30 minutes, and let stand for 4 hours or longer until the precipitate is decidedly crystalline. Filter into a weighed porcelain gooch, wash with  $(NH_4)_2HPO_4$ , ignite cautiously, until free from NH<sub>3</sub> and then at dull redness for 15 minutes. Weigh as  $Zn_2P_2O_7$ . 42.89 per cent. Zn.

Titration. The Zn may be titrated with 2.2 per cent. K<sub>4</sub>Fe(CN)<sub>6</sub>3H<sub>2</sub>O solution using 1 per cent. UO<sub>2</sub>(NO<sub>8</sub>)<sub>2</sub>6H<sub>2</sub>O solution, dropped on a wax coated tile, as an indicator. Standardize against 0.2 gram portions of ZnO dissolved in HCl. The conditions for this titration are discussed fully under Zn ores (pp. 60-62) and should be carefully studied before attempting a titration. Furthermore, it is well to try several titrations using a ZnCl<sub>2</sub> solution in order to acquire the technique of the process.

Determine the Zn in the filtrate from the Fe and Ni of the brass sample by one of the methods given above.

Analysis of	Sample No.	
		<del></del>
	,	
		•

i

#### **BRONZE.**

Sn only. Procedure No. 17.—The analysis of I gram of bronzes offers no great difference in procedure from the above except that the SnO<sub>2</sub> must be purified by fusion with a mixture of equal parts of S and Na<sub>2</sub>CO<sub>3</sub> = 6 times the weight of the SnO<sub>2</sub>. Dissolve the cooled melt in hot water. The polysulfids of the other metals are reduced with Na<sub>2</sub>SO<sub>3</sub> solution and the metals thus precipitating are filtered, washed and added to the main filtrate. The Sn remains in solution and may be precipitated as metal on a Cu plated Pt cathode. (Add 4 grams KCN to keep S in solution.) Or, better still, the foreign sulfids are ignited cautiously to oxids and weighed before dissolving and adding to the main solution. Their weight is deducted from the total weight of the impure SnO<sub>2</sub>.

Optional Rapid Procedure No. C.—The rapid analysis of Bronze is accomplished as follows: Dissolve I gram in aqua regia, boil off Cl, take up in 15 cubic centimeters HCl I-3, dilute to 150 cubic centimeters, add 2 grams NH<sub>2</sub>OH.HCl and electrolyze, keeping the cathode potential at —0.40 volt against a calomel electrode, using a centi or millivoltmeter. The Cu will deposit in about 40 minutes using rapid stirring and gauze electrodes and the end point is clearly seen by the increase in cathode potential. Whereupon the current is decreased and finally shut off.

Weigh the washed, and dried with  $C_2H_5OH$  and  $(C_2H_5)_2O$ , cathode. Difference between the initial and this weighing  $\times$  100 = per cent. Cu. Replace in the solution and after adding 3 grams more NH<sub>2</sub>OH.HCl electrolyze (40 minutes) until no Zn is found with H<sub>2</sub>S. In replacing the Cu plated cathode for the Sn deposit it should not set quite so deep in the solution as before. This prevents the Sn deposit coming directly in contact with the Pt and alloying with it.

In practice, most bronzes contain Zn (sometimes Pb) in addition to Cu and Sn, but the details of the rapid separation of these three or four elements from HCl solution have not yet been published.

Phosphor bronzes must be dissolved in HNO<sub>3</sub>, and the SnO<sub>2</sub> + impurities weighed. The Sn or impurities determined as under bronze, and the P determined as follows: Fuse the SnO<sub>2</sub> from 3-10-gram sample of the bronze with KCN to form potassium phosphate, boil under hood, rid of Sn with H<sub>2</sub>S, filter concentrate filtrate to 50 cubic centimeters, boil with Br water and precipitate the P with either

Analysis of	Sample No.	•	
		-	
	•		
	•		
			•

molybdate mixture from a HNO<sub>3</sub> solution or better with 10 cubic centimeters magnesia mixture\* per each 0.1 gram P<sub>2</sub>O<sub>5</sub> precipitated = 27.87 per cent. P.

Mn may be determined in bronze as follows: In absence of P dissolve bronze in HNO<sub>3</sub> filter, precipitate Pb as PbSO<sub>4</sub> details above, and electrolyze the filtrate. The Mn will deposit on the anode. Dissolve it off with HNO<sub>3</sub> and either perhydrol solution or a reducing agent, and determine the Mn as per the bismuthate method from the boiled HNO<sub>3</sub> solution. For other details and other bronzes, the larger books on the subject should be consulted.

German silver may be analyzed as follows: I gram is dissolved in HNO<sub>3</sub>, boiled, Cu and Mn precipitated by electrolysis, and Mn, dissolved from the anode, may be determined by the bismuthate method. The Ni in the electrolysis filtrate

is next determined as glyoxime, keeping the Fe in solution with HTr. Or the Fe may be separated by double precipitation with NH<sub>4</sub>OH, and the Ni (and Co) determined by electrolysis from the ammoniacal solution. These are separated by dimethylglyoxime followed by nitroso-beta-naphthol. The Fe is determined (after decomposing the HTr if used, with boiling HNO<sub>3</sub> or (NH<sub>2</sub>)<sub>2</sub>S<sub>2</sub>O<sub>8</sub> by NH<sub>4</sub>OH before or after the Ni (and Co if present). Finally the Zn is determined from the evaporated filtrate, by one of the methods described elsewhere.

Monel metal may be analyzed as follows: 0.5-gram dissolved in HNO<sub>3</sub> the Cu determined by electrolysis, the Fe precipitated by NH<sub>4</sub>OH, and the Ni determined by electrolysis in the filtrate. Or Cu may be separated with H<sub>2</sub>S, and determined by roasting to oxid, Ni determined as dimethylglyoxime, and the Fe determined in the filtrate after destroying the tartaric acid with boiling HNO<sub>3</sub> or  $(NH_4)_2S_2O_8$ . Other methods are described in Johnson.

When alloys contain large amounts of Ni, samples weighing only a few tenths of a gram should be taken for the dimethylglyoxime determination. If this precipitate weighs much more than 0.1-gram or if it gets cold, it becomes hard to handle. It is 20.31 per cent. Ni. Conditions for determination are given on p. 20.

## BEARING METALS."

The following procedure is recommended by Demorest. Zn must be determined in a separate sample after precipitation of the other elements with H<sub>2</sub>S.

Procedure No. 18.—Dissolve I gram of sample in 20 cubic centimeters H<sub>2</sub>SO<sub>4</sub> specific gravity 1.84 in a ½-liter Erlenmeyer, cool carefully, dilute to 70 cubic centimeters, reboil a few minutes and filter (preferably through a weighed alundum cone). The precipitate should be pure white, and the general conditions

<sup>\*</sup> Magnesia mixture. Eleven parts MgCl<sub>2</sub>6H<sub>2</sub>O, 28 parts NH<sub>4</sub>Cl both in 130 cubic centimeters H<sub>2</sub>O. To this add 7 grams NH<sub>4</sub>OH. Settle over night and decant or filter for use.

Analysis of	Sample No.		
-			
		•	
1			l .

under BaSO<sub>4</sub>, page 6, noted. The porcelain or alundum crucible containing the PbSO<sub>4</sub> should be suspended inside another. (Two small ring stand rings, the upper wound with bits of Pt wire where the alundum or porcelain touches are suitable.) Ignite the crucibles, the inner one becoming a dull red, continuing the ignition for about 25 minutes. PbSO<sub>4</sub> = 68.31 per cent. Pb. Or the Pb from PbSO<sub>4</sub> may be determined electrolytically as under Pb, p. 64.

Sb. Make up the filtrate to 120-140 cubic centimeters heat to about 80° C. and titrate about 0.5 cubic centimeter beyond the end point with 0.569 per cent. KMnO<sub>4</sub>. Titrate back with equivalent FeSO<sub>4</sub> solution, to the usual pink end point. (The last of the Sb is hard to oxidize, hence the excess KMnO<sub>4</sub>—mass action—is necessary. The excess is determined with FeSO<sub>4</sub> solution.) The KMnO<sub>4</sub> should be standardized against a warm dilute H<sub>2</sub>SO<sub>4</sub> solution of pure Sb, using equivalent FeSO<sub>4</sub> solution to titrate back with.

Cu. Cu is determined in the same solution after the KMnO<sub>4</sub> has been added, as follows: Add 3 grams HTr to keep Sb from precipitating, make slightly alkaline with NH<sub>4</sub>OH, add 2 cubic centimeters H<sub>2</sub>SO<sub>4</sub>, and heat to boiling. Reduce Cu<sup>++</sup> to Cu<sup>+</sup> with 2 grams powdered Na<sub>2</sub>SO<sub>3</sub>, and when the latter has dissolved, precipitate the Cu with 10 cubic centimeters of 10 per cent. KCNS solution. Filter, wash well with water containing a small amount of KCNS per liter, and ignite the Cu salt to CuO in a weighed porcelain crucible.

Other methods of Cu determination. Dissolve the CuCNS in HNO<sub>2</sub> and titrate with I solution. Dissolve the CuCNS in NaOH solution and titrate HCNS with KMnO<sub>4</sub> solution. Dissolve the ignited CuO and determine the Cu by electrolysis or by colorimetric methods.

Sn. The Sn content of the alloy may be determined as follows: Proceed as through Sb on another sample (the same sample if it contains no Cu), add 50 cubic centimeters HCl and transfer the solution to a 500 cubic centimeter flask, making up to about 200 cubic centimeters. This flask should be provided with a 4-hole rubber stopper. Illuminating gas or CO<sub>2</sub> bubbles through the solution and out through holes Nos. 1 and 3. The elongated tip of the burette passes through No. 2, and a rod of pure Fe, with the end bent into a coil passes through No. 4.

Analysis of	Sample No.		
		•	
	• •		
	•		

.

Set up the flask above a Bunsen. Boil while passing gas through the solution—in which the Fe coil dips—until the solution becomes colorless. Slip up the Fe ring from the solution, and titrate the reduced Sn with N/10 iodin solution using a little starch solution as an indicator. The iodin solution should be standardized against 0.2-gram Sn, to which 0.15-gram Sb has been added. The procedure is then the same as for the sample.

Iodine solution. 6.5 grams or 13 grams iodine (resublimed from KI) dissolved in 50 cubic centimeters H<sub>2</sub>O containing 12 or 20 grams KI diluted to 1 liter. Corresponding Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub> solution contains respectively 12.5 or 25 grams of thiosulfate.

A permanent starch solution may be made up by stirring 5 grams of starch into 200 cubic centimeters of cold water, heating to boil with constant stirring until the starch is thoroughly dissolved. Dilute with cold water to one liter and add 10 grams of crystallized ZnCl<sub>2</sub>. Settle over night and decant for use.

Starch solution. Add 1-gram starch plus 5 grams H<sub>2</sub>O after mixing thoroughly to 1 liter boiling water. Let stand, and use 2-3 cubic centimeters of the decanted solution.

Sn may be determined in tin plate by placing a definite fractional part of the total area of the plate in a tumbler  $\frac{2}{3}$  full of 40 per cent. Pb (Ac)<sub>2</sub> solution. The sample is dried and weighed before immersion, and after immersion the Pb carefully washed off. The final weight of the dried sample deducted from the initial  $\div$  the fraction of the total gives the weight of Sn on the sheet.

Sn may be separated from Sb (and Cu if present) by adding  $H_eS$  to a boiling oxalate solution of the metals for 45 minutes. Sn only remains in solution, and may be deposited on a copper plated cathode. 10 grams  $H_2C_2O_4$ , 10 grams  $(NH_4)$   $C_2O_4$  in 200 cubic centimeters  $H_2O$  form the solution.

Analysis of	Sample No.		
	-		
			•
·			
	<del></del>	·	

.

## MATTE AND ORES.

## MATTE.

**8**, see p. 8.

Pb, if present, will precipitate as sulphate on boiling off the acids and gases (see below) with H<sub>2</sub>SO<sub>4</sub>. It is determined according to the methods on pp. 38 and 62.

Cu. 0.5-1 gram of the matte is dissolved in 15 cubic centimeters HNO<sub>3</sub> specific gravity 1.42, pulverized KClO<sub>3</sub> added until all the sulfur is oxidized (to prevent inclusion of metals), and the solution boiled to beginning of fumes with 4.5 cubic centimeters H<sub>2</sub>SO<sub>4</sub>. Pb if present is filtered off from the diluted solution at this point and determined as per the preceding paragraph. The dilute Pb free solution is now electrolyzed for Cu after adding a few cubic centimeters HNO<sub>3</sub>.\* A few cubic centimeters H<sub>2</sub>O<sub>2</sub> solution are added toward the end of the electrolysis, as noted under brass.

Copper may also be determined by adding strips of Al to the Pb free solution after almost neutralizing (litmus) with NH<sub>4</sub>OH, and determining the copper as in steel, p. 22.

Fe, Zn. Iron is determined in the solution free from Cu by adding Al strips (or Zn if the latter (i. e., Zn) is to be determined) in a separate sample to reduce the Fe to the ferrous state. The iron content is then determined by KMnO<sub>4</sub> as on p. 14 ff. A blank must be run on the Zn and Al used as they both contain Fe; done by dissolving 10 grams in HCl and either titrating directly or else separating the Fe by Br and determining the Fe by titration or weighing as oxide.

Gold and silver are determined by fire methods. The determination of the lesser constituents is given in Classen and Price and Meade.

<sup>\*</sup>CuH<sub>2</sub>, copper hydrid, is formed momentarily on the cathode during electrolysis. HNO<sub>2</sub> decomposes (oxidizes) this more readily than H<sub>2</sub>SO<sub>4</sub>, hence the presence of a little HNO<sub>2</sub> is necessary in Cu analyses. Much HNO<sub>2</sub> prevents the precipitation of all the Cu when much Fe is present. See also p. 44, Sand ff.

Some authors suggest the direct electrolysis of a sulfate solution of matte, depending on the slight solubility of PbSO<sub>4</sub> to give finally all the Pb as PbO<sub>2</sub> on the anode.

Analysis of Sample No. **Zn.** Zinc may be determined as follows, using Al to reduce the iron, or using a separate sample.

Boil down the solution containing all some or none of the following (Fe, Al, Mn, Cu) and Zn to 10 cubic centimeters with HNO<sub>3</sub>, specific gravity, 1.42, add KClO<sub>3</sub> as above, and evaporate just to dryness. Cool, dissolve in concentrated HNO<sub>3</sub> and filter off Mn, on asbestos if present. To the filtrate, add 10 grams solid NH<sub>4</sub>Cl 15 cubic centimeters NH<sub>4</sub>OH, and 10 cubic centimeters Br water, boil, filter and wash with warm dilute NH<sub>4</sub>Cl solution. Filtrate contains (Cu) and Zn. Rid of the former, if present by passing H<sub>2</sub>S to the warm diluted slightly acid solution, or electrolyze. Rid the filtrate of H<sub>2</sub>S if used by boiling with HCl, and titrate the Zn with K<sub>4</sub>Fe(CN)<sub>6</sub> solution as on p. 42 or 60.

Analysis of	Sample No.		-3-
•			
		<del></del>	
	<del></del>		
			·
•			

,

## ORES.

## COPPER ORES.

There are six methods in use for the determination of Cu in ores; the electrolytic, iodid, sulphocyanid-permanganate, iodate, colorimetric and cyanid arranged in the probable order of greatest accuracy.

The first of these is the most acceptable as it involves only solution and two weighings. It has the disadvantage that it cannot be used with ores containing As, Sb or Bi. They may first be separated, p. 48, and the Cu determined electrolytically.

The iodid is the all-around method. It is not quite as accurate as the electrolytic, and is not used on refined copper or on pure ores in large laboratories.

The sulphocyanid allows of titration with KMnO<sub>4</sub> but is not as well established as the others. The iodate method is used by a few chemists.

The cyanid method is used for more approximate determinations.

For example, at one of the western plants the electrolytic method is used in the absence of As, Sb, Bi (or Pb). If separations are necessary, the iodid method is employed. For approximate results, the cyanid method is used, or for less than 0.5 per cent. Cu the colorimetric method is used. This last is used chiefly for slags and low-grade ores.

Procedure No. 19.—To determine Cu by the iodid method, I gram of ore is dissolved in ¼-liter Erlenmeyer, in HNO<sub>3</sub>, a few cubic centimeters of HCl added to rid of Ag, and the solution boiled with H<sub>2</sub>SO<sub>4</sub> to white fumes. Cool, add 30 cubic centimeters water and heat until Fe<sub>2</sub>(SO<sub>4</sub>)<sub>8</sub> dissolves. Precipitate the Cu with either Al or Zn (followed by a little H<sub>2</sub>S) dissolve all the Zn or Al, filter and wash Cu and insoluble on paper. Wash these back into the beaker, using but little water, and dissolve in 10 cubic centimeters HNO<sub>3</sub> 1-2. Boil the solution, pour through the filter, wash with a little Br water to dissolve any CuS caught mechanically, and boil off the Br 5-10 minutes. It has oxidized any As or Sb so that they will not interfere with the titration. Add alkali to first traces of blue color, and distinctly acidify with HAc., forming Cu (Ac)<sub>2</sub>. To this add 3 grams

	Analysis of	Sample No.		
				•
•				
				·
				<del></del>
				ı
			<del></del>	
			<del></del>	

.

KI which has been dissolved in a few cubic centimeters water. This forms CuI, liberating one equivalent of iodin for every equivalent of Cu present. This I is titrated at once with 1.96 per cent. Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub> solution standardized against a known ore, using starch solution, p. 52, as an indicator.

The procedure in the cyanid titration is the same up to the adding of alkali. Faint blue is produced with NH<sub>4</sub>OH, and the color almost discharged (giving a pink tint) with ore standardized 2.1 per cent. KCN solution.

A moisture test should be made on the ore in porcelain at 100° C.

### ZINC IN ORES."

The titration of Zn involves the non-academic procedure of titrating ¾ of the solution, until the end point is passed; then adding most of the remaining ¼ of the solution titrating to about the end point, and finally rinsing out the auxiliary beaker and getting the true end point. The indicator ordinarily used is a row of drops of I per cent. UO<sub>2</sub>(NO<sub>8</sub>)<sub>2</sub>6H<sub>2</sub>O solution on a paraffin coated porcelain tile. Breyer has suggested the use of the internal "blue color" in the solution which changes to green at the end point, when Fe<sup>II</sup> is present in the Zn solution and <sup>III</sup>Fe(CN)<sub>8</sub> in the ferrocyanid.

The Breyer method given below is for the analysis of complex ores. If Cd is not present, only one H<sub>2</sub>S treatment is necessary; if the substance is insoluble e. g., zinc spinel, it must be fused with Na<sub>2</sub>CO<sub>3</sub> and borax to dissolve it. If pure zinc from brass is to be titrated the requisite amount of NH<sub>4</sub>OH and HCl is added to X, p. 62, the 0.4 mg. Fe<sup>11</sup> added and the procedure thence followed:

Procedure No. 20.—0.5-1 gram of ore is weighed into a tall covered 250 cubic centimeter Jena beaker, moistened with water and dissolved in 10 cubic centimeters HCl 1.20. Boil the solution, add carefully 10 cubic centimeters H<sub>2</sub>SO<sub>4</sub> 1-1 and boil to fumes. Dilute to 150 cubic centimeters and add 1 gram Al powder. Boil to colorless solution. Filter, pass in H<sub>2</sub>S. SiO<sub>2</sub> and H<sub>2</sub>S metals are removed, solution contains Al, Fe, Cd, Zn (Cd is rid of by H<sub>2</sub>S or by electrolysis). Procedure; add H<sub>2</sub>S to solution, a few drops only of NH<sub>4</sub>OH to start precipitates and filter on a paper washed with excess of (NH<sub>4</sub>)<sub>2</sub>S<sub>2</sub>. Rid solution of H<sub>2</sub>S by boiling.

Analysis of Samp		
	l	
	·	

Almost neutralize with KOH, finish neutralization with Na<sub>2</sub>CO<sub>3</sub> (litmus), add a few cubic centimeters H<sub>2</sub>SO<sub>4</sub> 1-5 and precipitate Zn with H<sub>2</sub>S. (Zn may be determined from here by roasting and weighing as oxid.) Preferably wash the precipitate on a filter paper, punch a hole in the latter, wash ZnS into a beaker, dissolve in 30 cubic centimeters HCl 1-3, boil off H<sub>2</sub>S for 15 minutes, add 13 cubic centimeters NH<sub>4</sub>OH, specific gravity 0.90, neutralize this with 3 cubic centimeters HCl 1.20.

X (Brass) dilute to 200 cubic centimeters, and:

Titrate the Zn with 2.2 per cent. "aged"  $K_4$ Fe  $(CN)_6.3H_2O$  solution as per below. Use either the  $UO_2(NO_3)_26H_2O$  solution on tile with the "blue color" check or else the latter only as indicator. Standardize against freshly ignited 0.2-gram ZnO dissolved in 30 cubic centimeters HCl 1-3 + 13 cubic centimeters  $NH_4OH$ , specific gravity 0.90, and using the same indicator as for the ore.

Cause of the "blue color" indicator action. If the ferrocyanid solution is allowed to age, gases dissolved in the water will by oxidation cause the formation of <sup>III</sup>Fe (CN)<sub>6</sub>. This latter does not react with Zn, but with the Fe<sup>II</sup> in the Zn solution. When all Zn is precipitated, the concentration of Fe<sup>II</sup> CN<sub>6</sub> becomes so great that it causes the decomposition of the blue Fe<sub>3</sub> [Fe (CN)<sub>6</sub>]<sub>2</sub> forming the colorless Fe[Fe (CN)<sub>6</sub>]. This color change takes two or three drops at the end point, running from blue through green to colorless. The most satisfactory point to recognize is the change from blue to green. This change comes at about the same place (0.25 cubic centimeter before actually) as the brown color produced on the tile. Hence it may be used to help get the brown color point; or better still, as the actual indicator for both standard and sample.

As directed under brass, it is advisable to experiment with a 1 per cent. solution of ZnCl<sub>2</sub> until the technique of the divided titration and of the end point is mastered. If the ferrocyanid solution is new, 0.3-gram, ferricyanid may be added. If the solution titrated has no Fe, 3-4 milligrams Fe<sup>11</sup> may be added before titration as directed under brass above. Finally, the same amount of NH<sub>4</sub>OH and HCl should be present in all titrations, viz.: 10 cubic centimeters HCl, specific gravity 1.20, and 13 cubic centimeters NH<sub>4</sub>OH, specific gravity 0.90. This is important.

A moisture test at 100° C. porcelain should be made on the ore.

#### LEAD IN ORES."

Pb is most readily determined as PbO<sub>2</sub>, by the electrolysis of a nitrate solution of the ore, but in the presence of As, Bi, Mn, As (Se, Te, P<sub>2</sub>O<sub>5</sub>) the separa-

Analysis of Sample No.

tion is incomplete. Therefore Pb is separated by precipitation as SO4, and this precipitate is treated according to the molybdate, bichromate or by electrolysis method. Procedure is to dissolve 1-gram ore in 10 cubic centimeters HCl 1.2, rid of the latter with 15 cubic centimeters HNO<sub>3</sub> 1.42 by boiling to small bulk. Most of the excess HNO<sub>3</sub> should be boiled off. If the ore is pure, dilute to 150-200 cubic centimeters and electrolyze. If the ore is not pure, ie., if it contains one of the above mentioned impurities, Pb should be precipitated as PbSO4 after ridding of most of the HNO3 neutralizing the remainder, diluting to 50 cubic centimeters 3 grams HTr to keep Sb in solution (p. 5), and adding diluted  $H_2SO_4$ to it or vice versa (p. 6). Five cubic centimeters C2H5OH should be added to the precipitate and about 10 cubic centimeters added to the wash bottle when washing the PbSO<sub>4</sub>. The bulk of this should remain in the original beaker. If the PbSO<sub>4</sub> is to be electrolyzed it is dissolved in hot saturated NH<sub>4</sub>NO<sub>3</sub> solution, run through the paper onto the precipitate in the beaker, and the solution diluted before adjusting the gauze anode. The electrolysis procedure is the same as under brass and bearing metals.\*

For impure ores the molybdate or the bichromate titration is more often used.

Procedure No. 21.—Weigh 0.5-1 gram ore into a 1/4-liter Erlenmeyer, add 10 cubic centimeters concentrated HCl and 20 cubic centimeters water, and boil to rid of all H<sub>2</sub>S. Any insoluble residue other than SiO<sub>2</sub> is dissolved by the addition of 10 cubic centimeters IINO<sub>3</sub> and 7 cubic centimeters H<sub>2</sub>SO<sub>4</sub> added to precipitate the PbSO<sub>4</sub>. Boil to white fumes to rid of HNO<sub>3</sub>, stand to cool, dilute to 50 cubic centimeters and allow to stand until the Fe<sub>2</sub>(SO<sub>4</sub>)<sub>3</sub> dissolves. Cool filter, and wash, adding a few cubic centimeters H<sub>2</sub>SO<sub>4</sub> to the wash bottle.

Transfer paper with precipitate to original flask. Disintegrate the paper by boiling with 10 cubic centimeters HCl 1.20 total bulk about 70 cubic centimeters. Add 15 cubic centimeters more HCl, cool, add a little water if too dilute, and

<sup>\*</sup> The precipitation of  $PbO_2$  from pure solution is aided by the addition of a small amount of  $Cu(NO_3)_2$  solution. Cu++ discharges more readily than H+, hence the opposite ion  $PbO_2--$  will thus more readily precipitate.

Analysis of Sample No. make distinctly alkaline with NH<sub>4</sub>OH. Make distinctly acid with glacial HAc. Heat to dissolve completely the PbSO<sub>4</sub>, and dilute to 200 cubic centimeters with boiling water.

Divide the solution into two parts and titrate hot with 0.425 per cent. (NH<sub>4</sub>)<sub>2</sub>MoO<sub>4</sub> solution. The indicator used is 0.5 per cent. tannin solution dropped onto the tile. A yellow color develops at the end point. The same divided titration is used as with Zn; and it is advisable to study this, using a 0.15 per cent. solution of Pb (NO<sub>3</sub>)<sub>2</sub> before attempting the actual titrations. The molybdate solution is standardized against 0.2230-gram PbSO<sub>4</sub> treated as the ore beginning with the litmus and NH<sub>4</sub>OH treatment.

A moisture test in porcelain at 100° C. should be run on the ore.

### ALO, IN BAUXITE."

**Procedure No. 22.**—Fuse 5 grams of pulverized sample with Na<sub>2</sub>O<sub>2</sub>, dissolve the melt in H<sub>2</sub>O and re-fuse the insoluble in Na<sub>2</sub>O<sub>2</sub>, dissolve again and make up the combined filtered filtrates to 1 liter. Acidify an aliquot part usually 1/10 with HNO<sub>3</sub>, add a slight excess, NH<sub>4</sub>OH, redissolve, reprecipitate. Filter, ignite in Pt, treat with H<sub>2</sub>F<sub>2</sub> + H<sub>2</sub>SO<sub>4</sub>, ignite on the ring burner and weigh as Al<sub>2</sub>O<sub>3</sub> 53.03 per cent. Al.

A moisture test is advisable.

Analysis of Sample No.

		·	
•			
		•	
_			
	·		

# **BIBLIOGRAPHY.**

1. Sampling, Camp, Tr. 8th Int. Cong. App. Chem., vol. 1, p. 91, for pig iron.

Sampling of iron ores. Crowell and Murray, "The Iron Ores of Lake Superior," Chaps. VI and VIII.

Sampling of lead ores. Hofman, "Metallurgy of Lead."

Sampling of copper ores. E. D. Peters, "Practice of Copper Smelting," p. 8.

Cleaning bronze, etc. Treadwell, p. 236.

- 2. Modern electrolytic methods of rapid analysis.
  - a. Sand, J. Chem. Soc. (London), 91, pp. 373-93, 1572.
  - b. Hildebrand, "Some Applications of the Hydrogen Electrode in Analysis, etc."

    J. Am. Chem. Soc., 35, p. 847.
  - c. Forbes and Bartlett, "Electrometric Titration of Ferrous Salts," ibid., p. 1527.
  - d. Shoch, Tr. 8th Int. Cong. App. Chem., 21, p. 81. Shoch, see note 20.
- 3. Silica and silicon. In fire sands, etc., Blair, p. 293.

In ferros, Johnson.

In rocks, "The Analysis of Silicate and Carbonate Rocks," W. F. Hildebrand. (Reprinted as Bull. U. S. G. S., No. 422.)

4. Phosphorus. Separation using standard acetate method. Blair, p. 80.

In mineral waters, Treadwell, p. 447.

In bronzes, see Bronze, No. 20.

Emmerton, Tr. Am. Inst. Min. Eng., vol. 15, p. 93. (KMnO<sub>4</sub> method.)

On yellow precipitate, Baxter, Am. Chem. J., 38, p. 298.

By Mg precipitation, Drown and Shimer, Tr. Am. Inst. Min. Eng., vol. 10, p. 137.

5. Sulfur volumetrically, Blair, pp. 68 and 329; Johnson, p. 173. Tr. 8th Int. Cong. App. Chem., 1, p. 97, standard method.

In pyrites, new standard method, ibid., p. 33.

6. Carbon, Sources of error in carbon determinations. Circular No. 14, Bureau of Standards.

Combustion with volumetric estimation of CO<sub>2</sub>. Stahl u. Eisen, 33, p. 449.

Other less usual methods, see chapter in Lord and Demorest.

New colorimeter. White, J. Am. Chem. Soc., 34, p. 659.

Determination of carbon in steel and iron by the barium carbonate methods. Cain, J. Ind. Eng. Chem., vol. 6, p. 465.

Manganese. Details of bismuthate method and of standardization of KMnO. solutions.
 Reprints of Bureau of Standards, Nos. 186 and 182.

Discussion of and references to other methods. Circular, ibid., No. 26.

Cain's modification for alloy steels. J. Ind. Eng. Chem., 5, p. 647.

Details of other methods. Low, p. 162 ff.

Standard basic acetate method. Blair, p. 108.

8. Fe titration. Reprint No. 182, Bureau of Standards.

Am. Chem. J., 33, pp. 500-34, 109.

See also note 2, Hildebrand and Forbes.

Procedure in detail for ores. Low, Demorest, Blair.

Diphenylcarbazid inside indicator for K<sub>2</sub>Cr<sub>2</sub>O<sub>7</sub> titration. J. Am. Chem. Soc., 35, p. 156.

The permanganate determination of iron in the presence of chlorids. J. Am. Chem. Soc., 36, p. 1429.

9. Nickel. Commercial nickel assay. Ann. Chem. Anal., 18, p. 385.

Brunck, Stahl u. Eisen, 28, p. 331.

Prettner, Chem. Ztg., 33, p. 396.

Circular, Bureau of Standards, No. 14, p. 13.

By electrolysis, references. E. F. Smith, Electro-Analysis, pp. 130, 264.

Cyanid titration of glyoxime. Blair, p. 333.

Cobalt separation. Knorre u. Illinski Ber., 1885, p. 699. As nitroso beta naphthol.

Details of Ni and Co electro analysis. J. Am. Chem. Soc., 33, p. 493.

Details of ether separation, Blair, 7th ed., p. 191.

10. Copper. Circular, Bureau of Standards, No. 14, p. 12.

Price, J. Ind. Eng. Chem., vol. 6, p. 170.

11. Chromium. Circular, Bureau of Standards, No. 14, p. 13; ibid., Tech. Paper, No. 6.

Chromium in ores. Dept. of Mines, Canada. Report on Chrome Iron Ores Deposits, chap. 11, p. 97.

Rapid method for Cr and Cr Ni steels. Low, p. 82.

Lord and Demorest, p. 149 ff.

12. Vanadium. Circular, Bureau of Standards, Nos. 26 and 14.

By inclusion in yellow precipitate. Tech. Paper No. 8.

Separation of Fe and Cr by electrolysis. Reprint No. 161.

Determination of Mn, Va, and Ti in same sample. Wright, Chem. News, No. 108, p. 248.

Va and Ti in same sample. Wright, Chem. News, No. 108, p. 248.

13. Chromium and vanadium. Johnson, J. Am. Chem. Soc., 30, p. 1229. Amplified Crites, J. Ind. Eng. Chem., vol. 3, p. 575. (Persulphate may be used in place of the strong permanganate.)

Optional method. Demorest, J. Ind. Eng. Chem., 4, p. 895.

14. Molybdenum. Circular, Bureau of Standards, No. 14, p. 12.

Blair, p. 207; E. F. Smith, p. 164.

15. Titanium gravimetric procedures. Lord and Demorest, p. 157.

Tr. 8th Int. Cong. App. Chem., vol. 1, pp. 285 and 445.

(Colorimetric and volumetric.)

Hildebrand, "Analysis of Silicate and Carbonate Rocks," Bull. U. S. G. S., No. 422.

The volumetric determination of titanium and chromium by means of a modified Jones reductor. J. Am. Chem. Soc., 36, p. 1426.

16. Aluminum. Blair, pp. 198, 254.

Better methods. Circular, Bureau of Standards, No. 26; referring to Ether method, J. Ind. Eng. Chem., 2, p. 45.

As phenylhydrazine salt. Hess and Campbell, J. Am. Chem. Soc., vol. 121, p. 776; E. T. Allen, ibid., 25, p. 421.

Drown and McKenna, Tr. Am. Inst. Min. Eng., 20, p. 242.

E. F. Smith, Electro-Analysis, p. 255.

17. Limestones, etc. Hildebrand book lc, reprinted as Bull. U. S. G. S., No. 422.

Blair, p. 280 ff.

The analysis of blast furnace slag. J. Am. Chem. Soc., 19, p. 370.

18. Ca and Mg. Fox, J. Ind. Eng. Chem., 5, p. 910.

19. Mg. Hillebrand, Tr. 8th Int. Cong. App. Chem., vol. 1, p. 217.

20. Non-ferrous alloys, etc.

Cu, Sn alloy. J. Soc. Chem. Ind., 32, pp. 581-4.

Cu, Ni, Zn alloy. Chem. Ztg., 37, p. 1372.

German silver. Price and Meade, p. 197.

Bronze, complete analysis. Chem. Ztg., 37, p. 1464.

Brass electrodes. Met. Chem. Eng., Vol. 12, p. 87.

Separating metals by graded potential. See note 2a.

Bronze electrolysis from HCl solution. See note 2d.

Phosphorus in non-ferrous alloys. Treadwell p. 238.

Bismuth. Price and Meade, pp. 56, 176-8.

Volumetric determination of Hg. Ber. d. d. chem. Ges., 1913, p. 3853.

Volumetric determination of Sn. Eng. Min. Jr., 1914, p. 115.

Rapid analyses of brasses, of bronzes and of bearing metals. Cady and others, J. Ind. Eng. Chem., July, 1914.

Electro analysis of Pb, Sn, Cu, Sb from HCl solution. Schoch. Read, not yet published.

Details of Cu and Pb electro analyses. J. Ind. Eng. Chem., 2, pp. 195, 348.

Tin in ores. Chapter in Low and Demorest.

Spelter. Price and Meade, Lord and Demorest, Tr. 8th Int. Cong. App. Chem., 1, p. 183. See also note 21.

Evolution—Precipitation method for phosphorus. Price and Meade, p. 120.

#### 21. Methods for zinc.

Summary of methods. J. Ind. Eng. Chem., 4, p. 469.

Conditions for gravimetric phosphate determination. Dakin, Z. Anal. Chem., 39, p. 273.

Electrolytic method. Spear and others. Tr. 8th Int. Cong. App. Chem., 25, p. 97.

Frary, J. Am. Chem. Soc., Nov., 1907.

Ferrocyanid titration and critique of methods. Breyer, Tr. 8th Int. Cong. App. Chem., 25, p. 7.

For ores. Demorest, J. Ind. Eng. Chem., 5, p. 303.

Influence of Pb on ferrocyanid titration. Lehner, Tr. 8th Int. Cong. App. Chem., 1, p. 279.

"Blue Color" inside indicator. Breyer lc.

New method for Zn. Am. Inst. Metals Proceedings for Sept., 1914. (Precipitate as Zn Hg(CgS)4.)

22. Bearing metals. Demorest, J. Ind. Eng. Chem., 5, p. 842. See also Price and Meade, and Treadwell.

Volumetric determination of Sn. Eng. Min. J., 1914, p. 1151. Cady, note 20.

- 23. Pb in ores by KMnO<sub>4</sub>. Bollenbeck, Chem. Ztg., 33, p. 1142. Other methods in Low, Demorest and Betts.
- 24. Copper. Notes 10 and 20. The function of addition agents to the bath in the electrolytic deposition of copper. Tr. 8th Int. Cong. App. Chem., 21, p. 99.

Function of reducing agents. Transactions of the Faraday Society, IX, part 18, 2, pp. 11 and 186.

Exact assay of As and Sb in refined copper. Heath, J. Ind. Eng. Chem., 3, p. 75. In presence of much iron. Classen, pp. 25, 257 ff.

Accuracy attainable with the iodid method. J. Am. Chem. Soc., 34, pp. 422, 894. Cyanid method, details. See Low, "Technical Methods of Ore Analysis." Sulfocyanid permanganate method. Guess, J. Am. Chem. Soc., 24. Iodate method for Cu. Eng. Min. J. for 1914.

25. Analyses of Al in a booklet published by the Aluminum Co. of America.

-

## 89083896639



B89083896639A



ENGINEERING LIEFARY

30

